

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	09-623-5000	Revision: 1
Remarks		single use only

Specification

Product	Blender bags PE, 3.500ml, 380x510 mm, without filter, 25 pcs/bag, transparent, np pcr ready, sterile R
Dimensions	width: 380,00 mm (+/- 0,2 mm)
	length: 510,00 mm (+/- 0,2 mm)
	thickness: 70 µm
Colour	transparent
Nominal Capacity	3.500ml
Material Characteristics	polyethylene (PE)
Temperature range	-40°C - 80°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	free of detectable phthalate
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 500
	pcs per pallet: 14.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked

Approved by nerbe plus quality department