

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-223-0100	Revision: 1
Remarks		single use only

Specification

Product	Screw cap HDPE, for mailing container 44x114mm, with rim, natural
Colour	natural
Material Characteristics	high density polyethylene (HDPE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 250 pcs per carton: 250 pcs per pallet: 27.000
Minimum order qty.	250 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	