

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-243-1100	Revision: 1
Remarks		single use only

Specification

Product	Mailing container PE, 30x30x180 mm, with absorbing inlay and enclosed screw cap PE, for swabs, natural
Dimensions	inner diameter opening: 29,00 mm (+/- 0,2 mm)
	height: 85,00 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 6.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	