

# Technical Datasheet

Date: 28.07.17

Page 1 of 1



EN ISO 13485 certified

EN ISO certified		
<b>Article - No:</b>	<b>09-812-0001</b>	<b>Revision: 4</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	PS stick with small viscose tip, Ø2,5x150mm, np pcr ready, CE
<b>Classification</b>	CE
<b>Dimensions</b>	diameter: 2,50 mm (+/- 0,2 mm) length: 150,00 mm (+/- 0,2 mm) diameter tip : 4,00-5,50 mm (+/- 0,2 mm)
<b>Material Characteristics</b>	stick: polystyrene (PS) tip: viscose (VI)
<b>Product Characteristics</b>	biologically inert physiologically harmless free of detectable phthalate acc. to REACH regulation
<b>Purity Grade</b>	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
<b>Sterilisation Process</b>	N/A
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 100 pcs per carton: 9.000 pcs per pallet: 504.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable

Approved by nerbe plus quality department