

Technical Datasheet

Date: 28.07.17

Page 1 of 1

EN ISO certified		
Article - No:	09-819-0110	Revision: 2
Remarks		single use only

Specification

Product	Wooden applicator with small cotton tip, Ø2,2x230mm, np pcr ready, CE
Classification	CE
Dimensions	diameter: 2,20 mm (+/- 0,2 mm) length: 230,00 mm (+/- 0,2 mm) diameter tip: 5,00 - 5,50 mm (+/- 0,2 mm)
Material Characteristics	stick: wood (TI) tip: cotton (CO)
Product Characteristics	biologically inert physiologically harmless free of detectable phthalate acc. to REACH regulation
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	N/A
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 100 pcs per carton: 10.000 pcs per pallet: 360.000
Storage	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	