

Technical Datasheet

Date: 20.10.17

Page 1 of 1

EN ISO 13485 certified

EN ISO certified		
Article - No:	09-819-8002	Revision: 3
Remarks		single use only

Specification

Product	Wooden applicator with small cotton tip, Ø2,2x150mm, single peel-pack, np pcr ready, sterile EO, CE
Classification	CE
Dimensions	diameter stick: 2,20 mm (+/- 0,2 mm) length stick: 150,00 mm (+/- 0,2 mm)
Material Characteristics	stick: wood (TI) tip: cotton (CO)
Product Characteristics	biologically inert physiologically harmless free of detectable phthalate acc. to REACH regulation
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	gas sterilisation by ethyleneoxid
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100 pcs per carton: 1.000 pcs per pallet: 70.000
Storage	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	