

Technical Datasheet

Date: 11.11.16

Page 1 of 1

| EN ISO certified | | |
|----------------------|--|------------------------|
| Article - No: | 09-312-0000 | Revision: 2 |
| Remarks | max. vol.: approx. 11 l / working vol.: approx. 5 l | single use only |

Specification

| | |
|---------------------------------|---|
| Product | Autoclavable bag PP, 300x500 mm, thickness: 50µm, natural, autoclavable up to 134°C |
| Dimensions | width: 300,00 mm (+/- 8 mm) |
| | length: 500,00 mm (-10/+20 mm) |
| | thickness: 50 µm (+/- 10%) |
| Colour | natural |
| Material Characteristics | polypropylene (PP) |
| Autoclavability | autoclavable up to 134°C |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | extremely resilient double bottom seam |
| | density approx. 0,90 gr. |
| Packaging | pcs per cs: 1.000 |
| | pcs per carton: 1.000 |
| | pcs per pallet: 20.000 |
| Storage | in dry condition: 40% - 60% rel. humidity |
| | bags and cartons originally closed |
| | storage temperature +10°C to maximal +30°C |
| | pallets cannot be stacked |

Approved by nerbe plus quality department