

Technical Datasheet

Date: 11.06.18

Page 1 of 1

EN ISO certified		
Article - No:	09-332-0000	Revision: 3
Remarks	max. vol.: approx. 19 l / working vol.: approx. 6 l	single use only

Specification

Product	Autoclavable bag PP, 400x500 mm, thickness: 50µm, natural, autoclavable up to 134°C
Dimensions	width: 400,00 mm (+/- 8 mm)
	length: 500,00 mm (-10/+20 mm)
	thickness: 50 µm (+/- 10%)
Colour	natural
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	extremely resilient double bottom seam
	density approx. 0,90 gr.
Packaging	pcs per cs: 250
	pcs per carton: 250
	pcs per pallet: 25.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	