

Technical Datasheet

Date: 11.11.16

Page 1 of 1

EN ISO certified		
Article - No:	09-342-0700	Revision: 2
Remarks	max. vol.: approx. 32 l / working vol.: approx. 9,5 l	single use only

Specification

Product	Autoclavable bag PP, 500x600 mm, thickness: 50µm, natural, autoclavable up to 134°C
Dimensions	width: 500,00 mm (+/- 8 mm)
	length: 600,00 mm (-10/+20 mm)
	thickness: 50 µm (+/- 10%)
Colour	natural
Material Characteristics	polypropylene (PP)
Autoclavability	autoclavable up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	extremely resilient double bottom seam
	density approx. 0,90 gr.
Packaging	pcs per cs: 500
	pcs per carton: 500
	pcs per pallet: 12.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	