

Technical Datasheet

Date: 09.08.17

Page 1 of 2

EN ISO certified		
Article - No:	09-509-5000	Revision: 1
Remarks		single use only

Specification

Product	Test tube PP with fitted screw cap PE, 12ml, Ø16x100 mm, with wooden applicator with cotton ti tube transparent, cap blue, graduated, sterile R
Dimensions	outer diameter: 16,60 mm (+/- 0,2 mm)
	outer diameter (incl. cap): 19,60 mm (+/- 0,2 mm)
	length (cap and tube): 102,60 mm (+/- 0,2 mm)
	height cap : 14,40 mm (+/- 0,2 mm)
	cotton swab:
	diameter: 2,20 mm (+/- 0,2 mm)
	length: 85,00 mm (+/- 0,2 mm)
diameter tip: 5,00 mm (+/- 0,2 mm)	
Colour	tube transparent, cap blue
Graduation	grad.
Nominal Capacity	12ml
Material Characteristics	tube: polypropylene (PP)
	plug: polyethylene (PE)
Temperature range	-20°C - 60°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 1.000
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity