

Microscope Adapters for Micromanipulation Systems Calibre Scientific

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Cell Manipulation](#) / [Microscope Adapters for Micromanipulation Systems Calibre Scientific](#)

Product Data

Flexibility gives you the freedom you need. With the help of the different microscope adapters, you can fit your cell manipulation solution to your microscope of choice. Micromanipulation systems are compatible with most instruments of the major microscope manufacturers.

Product Information

Microscope adapters enable the simple installation of the TransferMan® 4r or InjectMan® 4 on all conventional inverted microscopes manufactured by Leica®, Nikon®, Olympus® and Zeiss®.

Item No.	type no.	description	Price
5192325007	Universal stand	for mounting Eppendorf micromanipulators on upright microscopes and stereo microscopes independent of the microscope tripod used	CHF 1'061.00

More product images





* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 23.05.2026 at 16:07 o'clock.

www.huberlab.ch