

Vial disposer Arnold / Marco, PP SCAT

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Disposal](#) / [Disposal containers](#) / [Vial disposer Arnold / Marco, PP SCAT](#)

Product Data

For the safe and quickly disposal of HPLC vials. Thanks to the safety funnel with an extra large neck diameter, your vials fall easily through the opening without jamming or clogging the opening. The hinged lid safely holds back solvent vapors and other contaminants.

- For the disposal of large quantities of sample vials
- Arnold: Connection for optional activated carbon filter available
- Marco: Low overall height, also suitable for safety storage cabinets

| Item No. | type no. | Volume | Dimensions | thread | diameter | Price |
|-------------|----------|--------|--------------------|--------|----------------|------------|
| 660.3073.90 | Arnold | 10 l | 195 x 195 x 570 mm | S90 | Funnel: 185 mm | CHF 401.35 |
| 660.3073.91 | Marco | 10 l | 180 x 180 x 450 mm | S90 | Funnel: 140 mm | CHF 315.00 |

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 26.05.2026 at 15:20 o'clock.

www.huberlab.ch