


Contact Dishes / Square Petri Dish

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Microbiology](#) / [Petri dishes](#) / [Contact Dishes](#) / [Square Petri Dish](#)

Product Data

Contact Dishes

Contact dishes are used in hygiene monitoring including detection of microorganisms and testing the effectiveness of cleaning and disinfection on flat surfaces by means of contact cultures. Particularly in the pharmaceutical and food industries, analysis of germ count and germ species on surfaces is of major importance to detect and identify possible sources of contamination.

- Sterilised by irradiation
- Vented and non-vented options
- Graduation and division into 10 mm squares to enable quick and easy analyses


CELLSTAR® OneWell Plate™

Filled with agar, the OneWell Plate™ can be used as an HTS petri dish for the cultivation and screening of bacteria. With an area of 95 cm², the OneWell Plate™ also fills the gap between the petri dishes Cat.-No. 7.664 102 (58 cm²) and Cat.-No. 7.639 102 (143 cm²) / 7.688 102 (144 cm²). Handling and the required incubator space are improved compared to a round cell culture dish. Notches on the left side of the plate and the lid ensure a secured lid position. Beside general bacteriological applications, the OneWell Plate™ can be used as a multipurpose liquid container or disposable for the denaturation, hybridisation and washing of membranes (Southern, Northern and Western Blot). For dot blotting and hybridisation the OneWell Plate™ can be used as a membrane holder and storage device.

Item No.	description	particularity	d x h	vents	packaging	Price
7.629 161	contact dish	graduated bottom	65 x 15 mm	-	30 packs x 20 pcs	CHF 235.58

7.688 102	petri dish	square	(120 x 120) x + 17 mm	24 packs x 10 pcs	CHF 238.79
7.688 161	petri dish	square	(120 x 120) x + 17 mm	24 packs x 10 pcs	CHF 257.40

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.05.2026 at 08:13 o'clock.