

# Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>11-931-6000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Flip top container PP, 10ml, Ø21x54 mm, flat bottom, natural, np pcr ready, sterile R, autoclavable up to 121°C
<b>Dimensions</b>	outer diameter (bottom): 21,00 mm (+/- 0,2 mm)
	outer diameter (top): 26,00 mm (+/- 0,2 mm)
	inner diameter (top): 18,40 mm (+/- 0,2 mm)
	outer diameter (lid): 33,60 mm (+/- 0,2 mm)
	height (without lid): 51,00 mm (+/- 0,2 mm)
	height (with lid): 54,00 mm (+/- 0,2 mm)
<b>Colour</b>	natural
<b>Nominal Capacity</b>	10ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-40°C - 70°C
<b>Autoclavability</b>	autoclavable up to 121°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Shelf - Life</b>	24 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 500
	pcs per carton: 1.000
	pcs per pallet: 30.000
<b>Minimum order qty.</b>	500 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked