

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-952-6000	Revision: 1
Remarks	filling level line at 33 ml	single use only

Specification

Product	Flip top container PP, 50ml, Ø29x85 mm, flat bottom, natural, np pcr ready, sterile R, autoclavable up to 121°C
Dimensions	outer diameter (bottom): 29,00 mm (+/- 0,2 mm)
	outer diameter (top): 31,20 mm (+/- 0,2 mm)
	inner diameter (top): 29,00 mm (+/- 0,2 mm)
	outer diameter (lid): 35,00 mm (+/- 0,2 mm)
	height (without lid): 83,00 mm (+/- 0,2 mm)
	height (with lid): 85,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	50ml
Material Characteristics	polypropylene (PP)
Temperature range	-40°C - 70°C
Autoclavability	autoclavable up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	24 months at controlled storage conditions
Packaging	pcs per bag: 600
	pcs per carton: 600
	pcs per pallet: 14.400
Minimum order qty.	600 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked