

Technical Datasheet

Date: 12.12.18

Page 1 of 1

EN ISO certified		
Article - No:	08-221-2005	Revision: 2
Remarks		single use only

Specification

Product	Universal specimen container with spoon PS, 20ml, Ø22x62 mm, with fitted pressure plug PE, container transparent, plug blue, CE/IVD
Classification	CE/IVD
Dimensions	diameter: 22,00 mm (+/- 0,2 mm) height: 62,00 mm (+/- 0,2 mm)
Colour	container transparent, plug blue
Nominal Capacity	20ml
Material Characteristics	container: polystyrene (PS) plug: polyethylene (PE)
Temperature range	-30°C – 55°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 250 pcs per carton: 1.000 pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable

Approved by nerbe plus quality department