

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	08-222-1000	Revision: 1
Remarks		single use only

Specification

Product	Universal specimen container PP, 15ml, Ø20x76mm, with enclosed screw cap PE, skirted, tube transparent, cap natural, graduated, max. RCF 4.000g, CE/IVD
Classification	CE/IVD
Dimensions	diameter: 20,00 mm (+/- 0,2 mm) height (without screw cap): 76,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	graduated
Nominal Capacity	15ml
Material Characteristics	container: polypropylene (PP) cap: polyethylene (PE)
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 1.000 pcs per carton: 1.000 pcs per pallet: 24.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	