

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	08-400-6000	Revision: 1
Remarks		single use only

Specification

Product	Sputum complete transport system, 24ml, content: inner container PS/PE 08-411-6000 and outer container PP/PE 08-422-1000, CE/IVD
Classification	CE/IVD
Dimensions	diameter: 79,00 mm (+/- 0,2 mm) height: 28,00 mm (+/- 0,2 mm)
Nominal Capacity	24ml
Material Characteristics	container: polystyrene (PS) plug: polyethylene (PE)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Sterilisation Process	radiation sterilization SAL 10 ⁻⁶ inner container and lid only
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 100 pcs per carton: 100 pcs per pallet: 4.800
Minimum order qty.	100 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	