

# Technical Datasheet

Date: 29.06.17

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>08-411-6000</b>	<b>Revision: 2</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Inner container PS for sputum complete system, 24ml, Ø70x23mm, with fitted screw cap PE, container transparent, cap natural, sterile R, CE/IVD
<b>Classification</b>	CE/IVD
<b>Dimensions</b>	diameter: 70,00 mm (+/- 0,2 mm) height: 23,00 mm (+/- 0,2 mm)
<b>Colour</b>	container transparent, cap natural
<b>Nominal Capacity</b>	24ml
<b>Material Characteristics</b>	container: polystyrene (PS) cap: polyethylene (PE)
<b>Temperature range</b>	-20°C - 60°C
<b>Product Characteristics</b>	biologically inert physiologically harmless no heavy metals according to EC/94/62
<b>Purity Grade</b>	do not contain any pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Packaging</b>	pcs per bag: 200 pcs per carton: 200 pcs per pallet: 9.600
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	