

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>08-502-4803</b>	<b>Revision: 1</b>
<b>Remarks</b>	<b>The survival rate of bacteria corresponds to the standard: DIN 58942-4 The filling amount of the liquid to be examined must be 10 ml, with a maximum tolerance of +/- 10%</b>	<b>single use only</b>

## Specification

<b>Product</b>	EXACTOBAC-U®, with fitted screw cap PE, 10ml, urine transport tube PP, Ø16x100 mm, with label, with lyophilized urine stabilizer, tube transparent, cap yellow, CE/IVD
<b>Classification</b>	CE/IVD
<b>Dimensions</b>	outer diameter: 16,00 mm (+/- 0,2 mm) height (withour screw cap): 100,00 mm (+/- 0,2 mm)
<b>Colour</b>	tube transparent, cap yellow
<b>Nominal Capacity</b>	10ml
<b>Material Characteristics</b>	tube: polypropylene (PP) cap: polyethylene (PE)
<b>Product Characteristics</b>	biologically inert physiologically harmless media: lyophilisate on boric acid sturdy, shock resistant high chemical resistance
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000 pcs per carton: 1.000 pcs per pallet: 24.000
<b>Minimum order qty.</b>	1.000 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
<b>Approved by nerbe plus quality department</b>	