

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-912-2100	Revision: 1
Remarks		single use only

Specification

Product	Hinged lid container PP, 80 ml, Ø45x80 mm, transparent, graduated at 80 ml, writing area, np pcr ready, autoclavable up to 121°C
Dimensions	outer diameter (without lid): 45,00 mm (+/- 0,2 mm)
	inner diameter: 43,00 mm (+/- 0,2 mm)
	height (without lid): 78,00 mm (+/- 0,2 mm)
	height (with lid): 80,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	graduated at 80 ml, writing area
Nominal Capacity	80ml
Material Characteristics	polypropylene (PP)
Temperature range	-20°C - 120°C
Autoclavability	autoclavable up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 350
	pcs per carton: 350
	pcs per pallet: 7.000
Minimum order qty.	350 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked

Approved by nerbe plus quality department