

Technical Datasheet

Date: 30.06.17

Page 1 of 2

		EN ISO certified
Article - No:	09-382-0000	Revision: 2
Remarks		single use only

Specification

Product	Sterilization double bag PP/LDPE, 37 liter, 600x900mm, thickness: 150 µm, outer bag yellow, inner bag transparent, autocl. up to 134°C
Dimensions	inner bag:
	width: 595,00 mm (+/- 5,00 mm)
	length: 870,00 mm (+/- 10,00 mm)
	thickness: 100,00 µm (+/- 13 %)
	outer bag:
	width: 600,00 mm (+/- 5,00 mm)
	length: 900,00 mm (+/- 10,00 mm)
	thickness: 50,00 µm (+/- 13 %)
Colour	outer bag yellow, inner bag transparent
Nominal Capacity	37 liter
Material Characteristics	inner bag: low density polyethylene (LDPE)
	outer bag: polypropylene (PP)
Autoclavability	autocl. up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	inner bag sealed with outer bag
	outer bag perforated upper part
	melting point inner bag: 87°C
Packaging	pcs per cs: 100
	pcs per carton: 100
	pcs per pallet: 2.400
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable