

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	09-385-0000	Revision: 1
Remarks		single use only

Specification

Product	Security closure PA, for sterilisation double bag 09-382-0000, natural, autoclavable
Colour	natural
Material Characteristics	polyamide (PA)
Autoclavability	autoclavable
Product Characteristics	biologically inert physiologically harmless
Packaging	pcs per bag: 1.000 pcs per carton: 1.000 pcs per pallet: 117.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	