

Technical Datasheet

Date: 14.02.17

Page 1 of 2

		EN ISO certified
Article - No:	09-511-5007	Revision: 1
Remarks	Complies with the requirement of MDD directive 93/42/EEC, class IIa: transient use as surgical invasive product is feasible.	single use only

Specification

Product	Culture swab, without transport medium, flexible PS stick with viscose tip, color code white, np pcr ready, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	tube: length: 145,00 mm (+/- 0,2 mm) outer diameter: 12,00 mm (+/- 0,2 mm) inner diameter: 10,70 mm (+/- 0,2 mm) plug: length: 20,00 mm (+/- 0,2 mm) outer diameter: 12,00 mm (+/- 0,2 mm) stick: length: 150,00 mm (+/- 0,2 mm) diameter: 2,50 mm (+/- 0,2 mm) diameter tip: 4,00 mm (+/- 0,2 mm) total length of tube and plug: 165,00 mm (+/- 0,2 mm)
Material Characteristics	tube: polypropylene (PP) plug: polyethylene (PE) swab: polysterene (PS) tip of swab: viscose (VI)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004 free of detectable phthalate acc. to REACH regulation free of latex
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	radiation sterilization SAI 10-6