

# Technical Datasheet

Date: 14.02.17

Page 1 of 2

EN ISO certified		
<b>Article - No:</b>	<b>09-556-8056</b>	<b>Revision: 1</b>
<b>Remarks</b>	<b>Complies with the requirement of MDD directive 93/42/EEC, class IIa: transient use as surgical invasive product is feasible.</b>	<b>single use only</b>

## Specification

<b>Product</b>	Transport swab with amies medium, aluminium stick with viscose tip, color code orange, single peel-pack, np pcr ready, sterile R, CE/IVD
<b>Classification</b>	CE/IVD
<b>Dimensions</b>	tube: length: 145,00 mm (+/- 0,2 mm) outer diameter: 12,00 mm (+/- 0,2 mm) inner diameter: 10,70 mm (+/- 0,2 mm) plug: length: 20,00 mm (+/- 0,2 mm) outer diameter: 12,00 mm (+/- 0,2 mm) stick: length: 150,00 mm (+/- 0,2 mm) diameter: 0,90 mm (+/- 0,2 mm) diameter tip: 1,50 mm (+/- 0,2 mm) peelpack: length: 38,00 mm (+/- 0,2 mm) width: 210,00 mm (+/- 0,2 mm) total length of tube and plug: 165,00 mm (+/- 0,2 mm)
<b>Material Characteristics</b>	tube: polypropylene (PP) cap: polyethylene (PE) swab: aluminium (AL) tip of swab: viscose (VI)
<b>Product Characteristics</b>	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004 free of detectable phthalate acc. to REACH regulation free of latex
<b>Purity Grade</b>	free of PCR inhibitors free of human DNA