

## Product Specification

# VIOLET PEARL

## Disposable Nitrile Examination Gloves

- Powder free
- Violet
- Chemical resistance protective gloves
- Latex-free
- Non-sterile



CE 0197



SAFETY MEETS QUALITY

### Material

100% Synthetic Nitrile Latex

### Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	2.5
Minor Visual defects	G-1	4.0

### Dimensions

Number	4801	4802	4803	4804	4805
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	76 ± 3 mm	84 ± 3 mm	94 ± 3 mm	105 ± 3 mm	113 ± 3 mm
Length	min. 240 mm	min. 240 mm	min. 240 mm	min. 240 mm	min. 240 mm

### Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.07 ± 0.02
Palm (in the middle of the palm)	0.07 ± 0.02
Fingertip (distance 15 mm to top of middle finger)	0.10 ± 0.02

### Physical Properties

Standard	Before Aging	After Aging
Tensile Strength (N)	min. 6.0	min. 6.0

### Powder

No powder lubricant added

### Protein content

No proteins contained

### Cuff

Beaded cuff

### Design

Ambidextrous

### Surface

Micro-roughened finger tips

### Color

Violet

### Quality assurance

EG-construction sample test TÜV Rheinland – CE 0197

PSA – EN 374 CAT III, EN 388, EN 420, 89/686/EEC

MPG – EN 455 Teil 1-4

EG-Standard 596/2009

EG-Standard 202/2014

EN ISO 9001:2008

EN ISO 13485:2012

### Packing Style

100 pcs per dispenser, 10 dispensers per carton

### Hint

The above information apply to the standard version.  
Subject to change.

HUBERLAB. AG  
Industriestrasse 123  
4147 Aesch

T 061 717 99 77  
F 061 711 93 42

www.huberlab.ch  
info@huberlab.ch

**HUBERLAB.**

committed to science