

Buffer solution traceable to SRM from NIST and PTB pH 4.00 (20°C) Certipur®,

[Product in eShop](#)

Home / [Brands](#) / [Merck](#) / Buffer solution traceable to SRM from NIST and PTB pH 4.00 (20°C) Certipur®,

Product Data

HS-No	3822 00 00
Lagerung	+15°C bis +25°C
LGK	10-13
WHC	nwg

Item No.	Content	Packaging	Price
1.09435.1000	1'000 milliliter	Kunststoffflasche	CHF 86.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 08.05.2026 at 23:52 o'clock.