

## WHEATON inserts



Product in eShop

Home / [Products](#) / WHEATON inserts

### Product Data

Item No.	Description	UNSPSC	ECLASS	Price
8.1191.11	WHEATON - 1.2 mL Clear Round Inserts, Blue Snap Caps with Pre-Cut FEP/Silicone Septa, 8 mm, Type I	41100000	32030000	CHF 282.30
8.1192.26	WHEATON - 350 uL Clear Round Inserts, Blue Snap Caps with Pre-Cut FEP/Silicone Septa, 8 mm, Type I	41100000	32030000	CHF 282.30
8.1231.12	WHEATON - 1.0 mL Clear Conical Inserts in Deep Well Microplate, Polypropylene, 1 Kit (Case of 1)	41100000	32030000	CHF 157.60
8.1232.27	WHEATON - 1.2 mL Clear Conical Inserts in Deep Well Microplate, Polypropylene, 1 Kit (Case of 1)	41100000	32030000	CHF 160.25
8.1233.33	WHEATON - 2.2 mL Clear Conical Inserts in Deep Well Microplate, Polypropylene, 1 Kit (Case of 1)	41100000	32030000	CHF 177.65
8.1234.48	WHEATON - 1.2 mL Blue Cap Insert PTFE/Silicone Septa Kit (Case of 1)	41100000	32030000	CHF 314.95

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 23.04.2026 at 09:56 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)