

WHEATON media bottle, clear



Product in eShop

[Home](#) / [Products](#) / WHEATON media bottle, clear

Product Data

| Item No. | Description | UNSPSC | ECLASS | Price |
|-----------|--|----------|----------|------------|
| 8.2194.15 | WHEATON - Media Bottle, 125mL Clear Non-Graduated Without 33-430 Screw Cap, Autoclavable (Case of | 41100000 | 32030000 | CHF 340.80 |
| 8.2194.17 | WHEATON - Media Bottle, 250mL Clear Non-Graduated Without 33-430 Screw Cap, Autoclavable (Case of | 41100000 | 32030000 | CHF 578.40 |
| 8.2194.19 | WHEATON - Media Bottle, 500mL Clear Non-Graduated Without 33-430 Screw Cap, Autoclavable (Case of | 41100000 | 32030000 | CHF 369.60 |
| 8.2194.35 | WHEATON - Media Bottle, 125mL Clear Graduated Without 33-430 Screw Cap, Autoclavable (Case of | 41100000 | 32030000 | CHF 376.80 |
| 8.2194.37 | WHEATON - Media Bottle, 250mL Clear Graduated Without 33-430 Screw Cap, Autoclavable (Case of | 41100000 | 32030000 | CHF 578.40 |

| | | | | |
|-----------|--|----------|----------|------------|
| 8.2194.39 | WHEATON - Media Bottle, 500mL Clear Graduated Without 33-430 Screw Cap, Autoclavable (Case of | 41100000 | 32030000 | CHF 442.80 |
| 8.2194.95 | WHEATON - Media Bottle, 125mL Clear, Non-Graduated With 33-430 Solid Black Phenolic Polyethylene | 41100000 | 32030000 | CHF 480.00 |
| 8.2194.97 | WHEATON - Media Bottle, 250mL Clear, Non-Graduated With 33-430 Solid Black Phenolic Polyethylene | 41100000 | 32030000 | CHF 705.60 |
| 8.2194.99 | WHEATON - Media Bottle, 500mL Clear, Non-Graduated With 33-430 Solid Black Phenolic Polyethylene | 41100000 | 32030000 | CHF 482.40 |
| 8.2195.75 | WHEATON - Media Bottle, 125mL Clear, Non-Graduated With 33-430 Solid Black Phenolic 14B Rubber Lined | 41100000 | 32030000 | CHF 427.20 |
| 8.2195.77 | WHEATON - Media Bottle, 250mL Clear, Non-Graduated With 33-430 Solid Black Phenolic 14B Rubber Lined | 41100000 | 32030000 | CHF 633.60 |
| 8.2195.79 | WHEATON - Media Bottle, 500mL Clear, Non-Graduated With 33-430 Solid Black Phenolic 14B Rubber Lined | 41100000 | 32030000 | CHF 519.60 |

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 18.05.2026 at 09:59 o'clock.