

WHEATON Graduated media bottle



Product in eShop

Home / [Products](#) / WHEATON Graduated media bottle

Product Data

Item No.	Description	UNSPSC	ECLASS	Price
8.2199.17	WHEATON - Clear Borosilicate Glass 250mL Graduated Lab 45 Media Bottle, without 45mm Screw Cap (Case	41100000	32030000	CHF 280.20
8.2199.19	WHEATON - Clear Borosilicate Glass 500mL Graduated Lab 45 Media Bottle, without 45mm Screw Cap (Case	41100000	32030000	CHF 349.20
8.2199.20	WHEATON - Clear Borosilicate Glass 1000mL raduated Lab 45 Media Bottle, without 45mm Screw Cap (Case	41100000	32030000	CHF 447.00
8.2199.21	WHEATON - Clear Borosilicate Glass 2000mL Graduated Lab 45 Media Bottle, without 45mm Screw	41100000	32030000	CHF 614.70

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 21:29 o'clock.