

WHEATON Media bottle Lab 45



Product in eShop

[Home](#) / [Products](#) / WHEATON Media bottle Lab 45

Product Data

Item No.	Description	UNSPSC	ECLASS	Price
8.2199.27	WHEATON - Media Bottle, 250mL Lab 45 Wide Mouth, Graduated With 45mm Polypropylene Screw Cap (Case	41100000	32030000	CHF 310.20
8.2199.29	WHEATON - Media Bottle, 500mL Lab 45 Wide Mouth, Graduated With 45mm Polypropylene Screw Cap (Case	41100000	32030000	CHF 378.60
8.2199.30	WHEATON - Media Bottle, 1000mL Lab 45 Wide Mouth, Graduated With 45mm Polypropylene Screw Cap (Case	41100000	32030000	CHF 396.60
8.2199.31	WHEATON - Media Bottle, 2000mL Lab 45 Wide Mouth, Graduated With 45mm Polypropylene Screw Cap (Case	41100000	32030000	CHF 523.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 12.05.2026 at 10:13 o'clock.

www.huberlab.ch