

WHEATON media bottles



Product in eShop

Home / [Products](#) / WHEATON media bottles

Product Data

Item No.	Description	UNSPSC	ECLASS	Price
8.2674.47	WHEATON - 1000 mL Clear Media Bottle without Cap, Clear Borosilicate Glass Type 1 (Case of 12)	41100000	32030000	CHF 242.40
8.2675.53	WHEATON - Clear Borosilicate Glass 1000mL Graduated Media Bottle, without 38- 430 Screw Cap	41100000	32030000	CHF 243.00
8.2676.68	WHEATON - Clear Borosilicate Glass 1000mL Non-Graduated Media Bottle, with 38-430 Solid	41100000	32030000	CHF 303.00
8.2677.74	WHEATON - Clear Borosilicate Glass 1000mL Non-Graduated Media Bottle, with 38-430 Solid	41100000	32030000	CHF 347.40

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 21:14 o'clock.