

WHEATON Glass weight



Product in eShop

Home / [Products](#) / WHEATON Glass weight

Product Data

Item No.	Description	UNSPSC	ECLASS	Price
8.3743.20	WHEATON - Small Glass Sinkers, for 2mm and 3mm ID Tubing Assembly (Case of 6)	41100000	32030000	CHF 249.60
8.3743.21	WHEATON - Large Glass Sinkers, for 6mm and 8mm ID Tubing Assembly (Case of 6)	41100000	32030000	CHF 252.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.05.2026 at 12:50 o'clock.