

## b.safe Canister GL45 Bohlender

[Product in eShop](#)

[Home](#) / [Brands](#) / [b.safe](#) / b.safe Canister GL45 Bohlender

### Product Data

Canister with thread GL45 (DIN45) for the collection of liquid waste.

Item No.	l x w x h mm	material	capacity	Price
94.R205.02	150 x 110 x 210	PE-HD	2,5 l	CHF 18.85
94.R205.05	195 x 150 x 270	PE-HD	5,0 l	CHF 18.65
94.R205.10	225 x 190 x 300	PE-HD	10,0 l	CHF 20.50

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 27.04.2026 at 10:30 o'clock.