

## b.safe Politainers Bohlener

[Product in eShop](#)

Home / [Brands](#) / [b.safe](#) / b.safe Politainers Bohlener

### Product Data

Foldable canister for the collection of liquid waste. Please order additional accessories like the carton box or the adaptor for Waste Caps GL 45 separately.

#### Applications

Folded to save space during delivery and storage, unfolds automatically during filling. UN-approval only if used with the appropriate carton box.

Item No.	material	thread	l x w x h	description	Price
94.R295.05	PE	38 x 3	outer: 181 x 181 x 195 mm	Politainer 5L	CHF 9.25
94.R295.10	PE	38 x 3	outer: 228 x 228 x 228 mm	Politainer 10L	CHF 11.15
94.R297.20	PE	S 60	outer: 285 x 285 x 295 mm	Politainer 20L	CHF 13.95
94.R299.05	-	-	inner: 181 x 181 x 181 mm	Carton 5L	CHF 7.45
94.R299.10	-	-	inner: 228 x 228 x 228 mm	Carton 10L	CHF 8.35
94.R299.20	-	-	inner: 285 x 285 x 285 mm	Carton 20L	CHF 9.25

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 23.04.2026 at 09:21 oclock.

[www.huberlab.ch](http://www.huberlab.ch)