

Lab Drying Rack for counter or mount on wall Heathrow Scientific



Product in eShop

Home / [Brands](#) / [Heathrow Scientific](#) / Lab Drying Rack for counter or mount on wall Heathrow Scientific

Product Data

Bringing color to a mundane task

Most Versatile Rack in the Market

- Channels designed for more useful configurations
- Use on counter or mount on wall
- Large and Small pegs for all glassware varieties
- One of a kind Drying Basket to collect small items

Designed to protect against breakage

- Peg angle allows for better hold - especially large, heavy pieces
- SureGrip Hook securely grips glassware into place
- Holds up to 35 kg

Rack made of durable, easy to clean polystyrene

Dimensions: 44.1 x 66 x 13.9 cm

Included:

- 72 Pegs (36 Small Green, 36 Large Blue)
- 2 Drying Baskets
- Drainage Tubing
- Mounting Guide and Hardware

Item No.

15.1207.91

Price

CHF 420.25

[More product images](#)



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 22.04.2026 at 18:55 o'clock.

www.huberlab.ch