

## Disposal bags, standard, PP

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microbiology](#) / [Bags](#) / Disposal bags, standard, PP

### Product Data

With the autoclavable bags made of PP, you can dispose of your contaminated labware. Thanks to the strong PP film (50 µm) and an extra broad seam at the bottom, the bags are particularly tear-resistant. The bag size 200 x 300 mm is delivered for both versions in a handy dispenser box.

Autoclavable up to 134 °C. For autoclaving, the bags must only be slightly closed. Not highly transparent.

Item No.	Volume	Width x Height	Épaisseur du film	Price
32.2162.00	2,5 l	200 x 300 mm	50 µm	CHF 22.50
32.2162.10	8 l	300 x 500 mm	50 µm	CHF 276.29
32.2162.20	60 l	600 x 800 mm	50 µm	CHF 84.08
32.2162.30	38 l	400 x 780 mm	50 µm	CHF 353.49
32.2162.40	110 l	700 x 1100 mm	50 µm	CHF 100.54
32.2162.50	110 l	700 x 1100 mm	50 µm	CHF 503.45

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 26.04.2026 at 11:23 o'clock.