

## Cell rack, PP LLG-Labware

[Product in eShop](#)

Home / [Labortops Q2 2026](#) / Cell rack, PP LLG-Labware

### Product Data

Ideal for storing up to 16 cuvettes with 10 mm width. Two-tier arrangement of the cuvette rows ensures clarity.

- Numbered positions
- Autoclavable at 121 °C
- Dimensions (W x D x H): 210 x 70 x 35 mm

Item No.	type no.	Price
600.6820.61	LLG-Cuvette rack	CHF 8.30 CHF 20.85

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 07.05.2026 at 13:18 o'clock.