

Cuvettes, LLG-Labware



Product in eShop

[Home](#) / [Labortops Q2 2026](#) / Cuvettes, LLG-Labware

Product Data

LLG-Cuvettes are most frequently used in analytical laboratories.

They have two clear windows, two side walls and are available as optical glass or quartz variants. The walls and the base are polished inside and frosted outside. The thickness of the windows and the base is 1.25 mm. The overall height of standard type cells is 45 mm or 48 mm for those fitted with a stopper. Each cuvette is supplied with a PTFE lid.

- Optical Glass: 335 to 2500 nm
- Quartz: 200 to 2700 nm

Item No.	Description	UNSPSC	ECLASS	Price
600.6815.13	LLG-Cuvettes 3.5ml, quartz glass path length 10 mm, with PTFE Lid pack of 2	41100000	32030000	CHF 57.50 CHF 90.55
600.6815.14	LLG-Cuvettes 17.5ml, quartz glass path length 50 mm, with PTFE Lid pack of 2	41100000	32030000	CHF 214.00 CHF 337.00
600.6815.15	LLG-Cuvettes 3.5ml, optical glass path length 10 mm, with PTFE Lid pack of 4	41100000	32030000	CHF 13.50 CHF 21.15
600.6815.16	LLG-Cuvettes 17.5ml, optical glass path length 50 mm, with PTFE Lid pack of 4	41100000	32030000	CHF 28.75 CHF 45.25

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 17.04.2026 at 13:12 o'clock.

www.huberlab.ch