

Volumetric flasks, borosilicate glass 3.3, class A, amber glass LLG- Labware

[Product in eShop](#)

[Home](#) / [Labortops Q1 2026](#) / Volumetric flasks, borosilicate glass 3.3, class A, amber glass LLG-Labware

Product Data

Conformity certified with batch ID.

- DIN EN ISO 1042
- With standard ground joint and PE stopper
- White ring mark and inscription

Item No.	for nominal volume	Ground Joint	tol. +/-	Price
600.6862.46	5 ml	NS10/19	±0.040 ml	CHF 25.50
600.6862.47	10 ml	NS10/19	±0.040 ml	CHF 25.50
600.6862.48	20 ml	NS10/19	±0.040 ml	CHF 27.30
600.6862.49	25 ml	NS10/19	±0.040 ml	CHF 27.30
600.6862.50	50 ml	NS14/23	±0.100 ml	CHF 35.90
600.6862.51	100 ml	NS14/23	±0.100 ml	CHF 40.70
600.6862.52	200 ml	NS14/23	±0.150 ml	CHF 54.00
600.6862.53	250 ml	NS14/23	±0.150 ml	CHF 54.00
600.6862.54	500 ml	NS19/26	±0.250 ml	CHF 90.60
600.6862.55	1000 ml	NS24/29	±0.400 ml	CHF 70.10
600.6862.56	2000 ml	NS29/32	±0.600 ml	CHF 92.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 21:38 o'clock.