

Weighing bottles with NS lid, Borosilicate glass 3.3 LLG-Labware

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Multiuse containers](#) / [Jugs & glass containers](#) / [Weighing bottles with NS lid, Borosilicate glass 3.3 LLG-Labware](#)

Product Data

With writing area.

Available in tall or low shape.

Item No.	capacity	form	diameter	height	Ground Joint	Price
600.3001.50	11 ml	low	40 mm	25 mm	NS34,5/12	CHF 14.95
600.3001.51	3 ml	low	25 mm	25 mm	NS19/12	CHF 13.90
600.3001.52	10 ml	low	35 mm	30 mm	NS29/12	CHF 14.25
600.3001.53	20 ml	low	50 mm	30 mm	NS45/12	CHF 14.85
600.3001.54	29 ml	low	60 mm	30 mm	NS50/12	CHF 17.45
600.3001.60	9 ml	tall	25 mm	40 mm	NS19/12	CHF 13.80
600.3001.61	19 ml	tall	30 mm	50 mm	NS24/12	CHF 14.35
600.3001.62	42 ml	tall	35 mm	70 mm	NS29/12	CHF 14.75
600.3001.63	66 ml	tall	40 mm	80 mm	NS34,5/12	CHF 15.25
600.3001.64	105 ml	tall	50 mm	80 mm	NS45/12	CHF 14.85

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 29.04.2026 at 13:11 o'clock.