

Agilent



Product in eShop

Home / [Products](#) / Agilent

Product Data

Item No.	Description	UNSPSC	ECLASS	Price
210.0001.76	FID jet, universal fit, 0.011 inch id	41100000	32030000	CHF 158.40
210.0015.16	Swagelok, 1.6 mm ID, PEEK, finger-tight fitting 2/pack. for 1/16in OD tubing	41100000	32030000	CHF 49.60
210.0179.73	Gas Clean carrier gas purifier	41100000	32030000	CHF 492.80
210.0254.00	Sample rack for 16 mm od tubes, 60 positions	41100000	32030000	CHF 88.00
210.0403.11	Agilent column nut for Varian/Bruker, brass, 2/pk, Similar to Bruker 394955100	41100000	32030000	CHF 59.20
210.0511.75	Seal wash PTFE for 1260, 1290, 1220, 1200, 1100 pumps, 1 /pk	41100000	32030000	CHF 136.00
210.1001.92	Vial septa, 20 mm, tan PTFE/white silicone, 1000/cs. Septum size: 20 mm (for 20 mm cap)	41100000	32030000	CHF 848.00
210.1005.43	Agilent vial kit, 2 mL, 9 mm, amber, 100/pk. Vial size: 12 x 32 mm (12 mm cap)	41100000	32030000	CHF 124.80
210.1207.00	FID collector insulator	41100000	32030000	CHF 52.80
210.2007.15	2 ml screw cap sample vials, clear, labelling field, 100 pcs.	41100000	32030000	CHF 48.00
210.2560.01	A portfolio reference standards for United States - State specific methods.	41100000	32030000	CHF 284.80

210.2600.20	Cartridge for Active Inlet Valve, 600 bar	41100000	32030000	CHF 926.40
210.4008.20	InfinityLab long-life deuterium DAD/MWD lamp, with RFID tag. For Agilent diode array	41100000	32030000	CHF 1'936.00
210.4121.68	Glass filter, solvent inlet, 20 µm, for Agilent 1200 Infinity Series	41100000	32030000	CHF 120.00
210.5600.61	Filament, high temperature, EI ion source	41100000	32030000	CHF 254.88
210.6285.62	Cartridge for Active Inlet Valve, 400 bar	41100000	32030000	CHF 513.60
210.6365.86	Piston, sapphire, 100 µL, for 1100, 1200, 1260, 1120, 1220 pump and autosampler	41100000	32030000	CHF 355.20
210.6365.89	Piston seal, PTFE	41100000	32030000	CHF 264.00
210.6747.28	Seal cap. For all manual purge valves	41100000	32030000	CHF 76.80
210.6759.66	InfinityLab Quick Turn LC fitting	41100000	32030000	CHF 254.40
210.8088.53	Ferrule, graphite, short, 0.5 mm id, for 0.1 to 0.32 mm columns, 10/pk	41100000	32030000	CHF 116.80
210.8112.10	Cap, crimp, silver aluminum, 11 mm, PTFE/red rubber septa, 100/pk. Cap size: 11 mm	41100000	32030000	CHF 40.00
210.8112.11	Cap, crimp, silver alu., 11 mm, PTFE/silicone/PTFE septa, 100/pk. Cap: 11 mm	41100000	32030000	CHF 54.40
210.8112.17	Cap, crimp, red aluminum, 11 mm, PTFE/red rubber septa, 100/pk. Cap size: 11 mm	41100000	32030000	CHF 46.40
210.8112.70	Vial insert, 250 µL, glass with polymer feet, 100/pk Insert size: 5.6 x 30 mm	41100000	32030000	CHF 227.20
210.8133.75	Vial, crimp top, clear, certified, 2 mL, 100/pk. Vial size: 12 x 32 mm (11 mm cap)	41100000	32030000	CHF 32.00
210.8133.76	Vial, crimp top, amber, write-on spot 2 mL, 100/pk. Vial size: 12 x 32 mm (11 mm cap)	41100000	32030000	CHF 38.40
210.8133.77	Vial insert, 400 µL, glass, flat bottom, 500/pk Insert size: 5.6 x 31 mm	41100000	32030000	CHF 147.20

210.8205.40	Cap, snap, clear, PTFE septa, 100/pk. Cap size: 11 mm	41100000	32030000	CHF 41.60
210.8205.44	Vial, snap top, 2 mL, clear, 100/pk. Vial size: 12 x 32 mm (11 mm cap)	41100000	32030000	CHF 41.60
210.8205.45	Vial, snap top, 2 mL, amber, write-on spot, 100/pk. Vial size: 12 x 32 mm (11 mm cap)	41100000	32030000	CHF 46.40
210.8205.50	Cap, snap, clear, clear PTFE/red silicone septa, 100/pk. Cap size: 11 mm	41100000	32030000	CHF 46.40
210.8205.67	Vial, crimp/snap top, polypropylene, 1 mL, 100/pk. Vial size: 12 x 32 mm (11 mm cap)	41100000	32030000	CHF 46.40
210.8207.14	Vial, screw top, clear, certified, 2 mL, 100/pk. Vial size: 12 x 32 mm (12 mm cap)	41100000	32030000	CHF 44.80
210.8207.16	Vial, screw top, amber, write-on spot 2 mL, 100/pk. Vial size: 12 x 32 mm (12 mm cap)	41100000	32030000	CHF 51.20
210.8207.17	Cap, screw, blue, PTFE/red silicone septa, 100/pk. Cap size: 12 mm	41100000	32030000	CHF 49.60
210.8207.31	Septa, for screw top caps, PTFE/red silicone, 100/pk. Septum size: 9 mm (for 12 mm cap)	41100000	32030000	CHF 32.00
210.8208.37	Vial, headspace, crimp, flat bottom, 20 mL, 100/pk. Vial: 22.75 x 75 mm (20 mm cap)	41100000	32030000	CHF 86.40
210.8208.65	Vial pack, screw top, pre-assembled, clear vials with write-on spot, blue caps,	41100000	32030000	CHF 136.00
210.8208.66	Vial pack, screw top, pre-assembled, clear vials with write-on spot, blue caps,	41100000	32030000	CHF 147.20
210.8227.07	Frit, PTFE, 5/pk	41100000	32030000	CHF 80.00
210.8234.59	Cap, snap, red, clear PTFE/red silicone septa, 100/pk. Cap size: 11 mm	41100000	32030000	CHF 48.00
210.8320.76	Cap, screw, blue, preslit PTFE/silicone, 100/pk. Cap size: 12 mm	41100000	32030000	CHF 88.00

210.8344.77	Cap, crimp, aluminum, PTFE/silicone septa, 20 mm, 100/pk. Cap: 20 mm	41100000	32030000	CHF 139.20
210.8344.78	Cap, crimp, headspace, with septa, 20 mm, silver aluminum cap, PTFE/silicone septa,	41100000	32030000	CHF 147.20
210.8346.47	Inlet liner, split, single taper, glass wool, deactivated, low pressure drop	41100000	32030000	CHF 97.60
210.8853.65	Inlet liner O-ring, non-stick fluorocarbon, 10/pk	41100000	32030000	CHF 35.20
210.8865.94	Vial, snap top, with fixed insert, amber, 300 µL insert volume, 100/pk. Vial: 12 x 32 mm	41100000	32030000	CHF 212.80
210.9040.06	Inlet liner, Ultra Inert, splitless, dimpled, 2 mm id, 5/pk	41100000	32030000	CHF 345.60
210.9158.51	AVF 45 platinum, 1 qt	41100000	32030000	CHF 169.60

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 27.04.2026 at 10:31 o'clock.