

Cylinder with knobbed lid Duran® DWK

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sample preparation](#) / [Glassware](#) / [Cylinder with knobbed lid Duran® DWK](#)

Product Data

- Manufactured from borosilicate glass 3.3
- With polished rim
- The lid can be opened easily while using the button
- Fired-on, highly durable white ceramic

Ideally, containers that can be closed are used to store substances and materials. In this way, substances and materials can be protected from unwanted effects.

Item No.	capacity	diameter	height	Price
8.5010.01	250 ml	80 mm	80 mm	CHF 853.00
8.5030.01	500 ml	100 mm	100 mm	CHF 1'074.50
8.5050.01	1000 ml	120 mm	120 mm	CHF 114.85
8.5100.01	2000 ml	150 mm	150 mm	CHF 157.50
8.5210.01	6000 ml	210 mm	210 mm	CHF 305.15
8.5320.01	12000 ml	260 mm	260 mm	CHF 507.75

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 22.04.2026 at 18:55 o'clock.