


Dissolved oxygen bottles, Winkler pattern Starke



Product in eShop

Home / [Labortops Q4 2025](#) / Dissolved oxygen bottles, Winkler pattern Starke

Product Data

Winkler pattern. Soda-lime glass.

For determination of dissolved oxygen in water. The measured capacity is specified to $\pm 0,01$ ml. With white labelling area. Solid, bevelled glass stopper, secured with an accessory spring clip. Each bottle is adjusted to the relevant stopper. Stoppers and bottles are therefore not interchangeable. Each bottle and its stopper is marked with a unique, matching identification number.

Item No.	volume	NS cone	Price
683.3040.38	100 to 150 ml	14/23	CHF 26.90
683.3040.48	250 to 300 ml	19/26	CHF 31.45

Accessories

Item No.	Description	Price
5.6055.00	Clamps for oxygen determination bottles acc. to Winkler Lenz	from CHF 9.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 29.04.2026 at 13:11 o'clock.