

Grippers for vaccine vials, PP Heathrow Scientific

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Laboratory tools](#) / [Tongs](#) / Grippers for vaccine vials, PP Heathrow Scientific

Product Data

For precise gripping and safe removal of crimped vaccine/serum vials from cold storage boxes. Set consisting of 5 different colours (red, purple, blue, green, and yellow).

- Ergonomic squeeze area
- For vaccine vials up to Ø 21 mm
- Autoclavable at 121 °C

Dimensions (WxLxH): 20 x 152 x 35 mm

Item No.	type no.	Price
15.1209.09	Vaccine vial gripper (1x red, 1x purple, 1x blue, 1x green, 1x yellow)	CHF 109.85

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 10.05.2026 at 03:25 o'clock.