

Staining jars, soda-lime-glass LLG-Labware

[Product in eShop](#)

Home / [Labortops_Q1_2026](#) / Staining jars, soda-lime-glass LLG-Labware

Product Data

Schifferdecker staining jar:

- For 20 microscope slides 76 x 26 mm, horizontal

Hellendahl staining jar:

- For 16 microscope slides 76 x 26 mm, vertical

Makro staining jar:

- For 20 microscope slides with 76 x 26 mm (please order rack separately)

Rack, stainless steel, for Makro staining jar:

- With long handle, for 20 microscope slides 76 x 26 mm, horizontal

Each staining jar is supplied with a glass lid.

| Item No. | Description | UNSPSC | ECLASS | Price |
|-------------|--|----------|----------|-----------|
| 600.1632.50 | LLG-Staining jars Schifferdecker for 20 slides of 76 x 26 mm | 41100000 | 32030000 | CHF 12.30 |
| 600.1632.51 | LLG-Staining jars Hellendahl for 16 slides of 76 x 26 mm with extension | 41100000 | 32030000 | CHF 13.90 |
| 600.1632.52 | LLG-Staining jar Macro for 20 slides of 76x26mm | 41100000 | 32030000 | CHF 29.90 |
| 600.1632.53 | LLG-Staining rack, stainless steel for 20 slides, for LLG-staining jar Macro | 41100000 | 32030000 | CHF 52.85 |

[More product images](#)



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 22.05.2026 at 02:52 o'clock.

www.huberlab.ch