

# Roller Stool SuperStep Werner Dorsch

[Product in eShop](#)

[Home](#) / [Labortops\\_Q1\\_2026](#) / Roller Stool SuperStep Werner Dorsch

## Product Data

The treads made of TTP are not glued or printed, but injected, meaning they are inseparably connected to the body. The non-slip material makes slipping practically impossible.

- Optimisation of the roller bar geometry
- Easy-care surface
- 10-year warranty on the body
- Repair option (roller replacement) guaranteed
- 100% recyclable
- Made in Germany
- Height without load approx. 42.8 cm
- Height with load approx. 41.5 cm
- Diameter at top approx. 29.2 cm
- At bottom approx. 42.4 cm
- Weight: 2.15 kg
- Max. loading capacity: 150 kg

Item No.	colour	Price
698.2123.00	white	CHF 86.70
698.2123.01	black	CHF 86.70
698.2123.02	red	CHF 86.70
698.2123.03	blue	CHF 75.20
698.2123.12	grey	CHF 86.70
698.2123.37	light grey	CHF 86.70

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 19.05.2026 at 15:12 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)