

Test tube racks Switch-Grid, POM Bel-Art

[Product in eShop](#)

Home / [Labortops_Q1_2026](#) / Test tube racks Switch-Grid, POM Bel-Art

Product Data

This split-frame rack holds up to two different tube sizes at the same time. The grids can be switched quickly. The test tube rack can be easily cleaned with common lab detergents, ethanol, or isopropyl alcohol.

- For test tubes 10 mm to 30 mm diameter
- With 5 colour-coded grid inserts
- Three-tier design
- Non-slip rubber feet
- Chemical resistant
- Frame and grids are autoclavable at 121°C
- Material: POM
- Dimensions: 318 x 67 x 111 mm

Item No.	type no.	Price
15.1874.50	Switch-Grid	CHF 85.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.05.2026 at 05:13 oclock.

www.huberlab.ch