

Disposal cartons for broken glass Heathrow Scientific



Product in eShop

Home / [Labortops Q2 2026](#) / Disposal cartons for broken glass Heathrow Scientific

Product Data

Easy-access collection containers make it convenient to safely dispose of non-biohazardous broken glass from your lab.

- Sturdy corrugated container lined with a clear, thick plastic bag to protect against leakage
- Warnings printed in both English and Spanish for added safety and convenience
- Reclose-A/Remove-A-Lid allows you to reclose the container
- A wide-mouth opening easily accommodates disposal of large items
- Inner flap seals the lid when the container is filled
- Containers ship and store flat and easily fold together for use

Item No.	type no.	Volume	Price
15.2000.60	Bench Top Model	5 l	CHF 145.00 CHF 290.00
15.2000.61	Floor Model	63 l	CHF 85.00 CHF 170.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 17.04.2026 at 13:12 o'clock.

www.huberlab.ch