



Stativstangen, Stahl verzinkt Bochem



Product in eShop

Home / [Labortops Q2 2026](#) / Stativstangen, Stahl verzinkt Bochem

Product Data

Item No.	diameter	length	Price
13.5160.81	12 mm	500 mm	CHF 13.55 CHF 19.35
13.5161.82	12 mm	600 mm	CHF 14.95 CHF 21.40
13.5162.83	12 mm	750 mm	CHF 16.30 CHF 23.30
13.5163.84	12 mm	1000 mm	CHF 18.60 CHF 26.50
13.5164.85	12 mm	1250 mm	CHF 20.95 CHF 29.95
13.5165.86	12 mm	1500 mm	CHF 23.85 CHF 34.05

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 19.05.2026 at 10:11 o'clock.